

ERGO COMFORT® GEL


Therapy Platforms®

Anti-Fatigue Footwear Insoles

Patent #7913423 & other patents pending



Ergo Comfort **GEL** Shoe Insole



Dual Performance Cushioning

Full length **Ergosert™ technology** layer & added **ErgoGel™** shock pads provide superior cushioning and shock absorption to fight fatigue.

Extra Deep Heal Cup

Improves the foot's ergonomic cushioning and provides additional stability.

Anti-Blister Top

Keeps feet cooler, drier, & less prone to fatigue and blisters.

Price USD: \$ 41.95

- **Fatigue Fighting Cushioning ★★★★★**
Ergosert™ technology is laboratory certified to have the same ergonomic properties as premium industrial anti-fatigue matting. **Ergosert™** material is lightweight and durable, it does not harden like PVC foam, does not breakdown like latex foam or take a compression set like EVA foam.
- **Extra Deep Heal Cup**
Helps cushion the foot naturally, utilizing the cushioning tissues in the heel of the foot.
- **Tapered Toe Design**
Provides better fit in all shoes.
- **Air Flow**
Vents molded into the shape of the arch help move warm air out and draw cooler air into footwear.
- **Anti-Blister Top Cover**
By controlling the amount the foot slides across the insole's surface, the anti-blister top cover helps prevent blisters and keeps feet cooler & drier.
- **Versatile Fit and Use**
Ergobuddy® insoles will fit boots, hiking, casual, and athletic shoes. They are performance insoles to be used all day for standing, walking, or running.

PVC & Latex Free Insoles

Fitting: Remove the existing insole and slide the new Ergo Comfort Insoles into your shoes. If your new insoles need to be trimmed, use a pen and trace the old insole's shape (fore-foot part only) on the bottom side of the new insole. Carefully cut the new Ergo Comfort insoles with scissors and place them in your shoes.

Ergobuddy®, **Ergosert™**, and **ErgoGel™** are trademarks of Johnson Technologies Corporation.

SATISFACTION GUARANTEED

This product is guaranteed for one year against defects. If you have questions or comments, please contact Johnson Technologies Corporation
2034 Pittway Drive
Nashville, Tenn 37207
1-800-229-1930

www.ergobuddy.com

Ergosert™ material is Made in the USA

Applying the science of **ergonomics** to footwear™

A Step Ahead!